

2,3-dihydroxy-N-methyl-3-phenyl-N-[(1Z)-2-phenylethenyl]propanamide

Chemical Properties


CAS No. :

Formula: C₁₈H₁₉NO₃

Molecular Weight: 297.35

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	2,3-dihydroxy-N-methyl-3-phenyl-N-[(1Z)-2-phenylethenyl]propanamide is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.363 mL	16.8152 mL	33.6304 mL
5 mM	0.6726 mL	3.363 mL	6.7261 mL
10 mM	0.3363 mL	1.6815 mL	3.363 mL
50 mM	0.0673 mL	0.3363 mL	0.6726 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481