


## Setosusin

## Chemical Properties

CAS No. :	182926-45-0
Formula:	C <sub>29</sub> H <sub>38</sub> O <sub>8</sub>
Molecular Weight:	514.615
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Setosusin is a meroterpenoid fungal metabolite originally isolated from <i>C. setosus</i> . It reduces neurotoxicity induced by amyloid- $\beta$ (A $\beta$ ) aggregates in PC12 cells (EC <sub>50</sub> = 112.6 $\mu$ M).
Targets(IC50)	Others,Endogenous Metabolite
In vivo	In vivo, setosusin (30 mg/kg) induces tremors in mice.

## Solubility Information

Solubility	DMSO: Soluble (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9432 mL	9.7159 mL	19.4318 mL
5 mM	0.3886 mL	1.9432 mL	3.8864 mL
10 mM	0.1943 mL	0.9716 mL	1.9432 mL
50 mM	0.0389 mL	0.1943 mL	0.3886 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481