


## 2'-Deoxy-5-nitrocytidine

## Chemical Properties

CAS No. :	69100-02-3
Formula:	C <sub>9</sub> H <sub>12</sub> N <sub>4</sub> O <sub>6</sub>
Molecular Weight:	272.21
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	2'-Deoxy-5-nitrocytidine is a DNA methyltransferase inhibitor used in cancer research [1].
Targets(IC50)	Others,DNA Methyltransferase

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.6736 mL	18.3682 mL	36.7363 mL
5 mM	0.7347 mL	3.6736 mL	7.3473 mL
10 mM	0.3674 mL	1.8368 mL	3.6736 mL
50 mM	0.0735 mL	0.3674 mL	0.7347 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481