

Kushenol F

Chemical Properties

CAS No. :

Formula: C₂₅H₂₈O₆

Molecular Weight: 424.493

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Kushenol F is a useful organic compound for research related to life sciences and the catalog number is T124984.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3558 mL	11.7788 mL	23.5577 mL
5 mM	0.4712 mL	2.3558 mL	4.7115 mL
10 mM	0.2356 mL	1.1779 mL	2.3558 mL
50 mM	0.0471 mL	0.2356 mL	0.4712 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481