


## 3'-Deoxy-N6-isopentenyl adenosine

## Chemical Properties

CAS No. :	1237496-62-6
Formula:	C <sub>15</sub> H <sub>21</sub> N <sub>5</sub> O <sub>3</sub>
Molecular Weight:	319.36
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	3'-Deoxy-N6-isopentenyl adenosine is a valuable organic compound for life sciences research (Catalog No.: TNU0457; CAS No.: 1237496-62-6).
Targets(IC50)	Nucleoside Antimetabolite/Analog,Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.1313 mL	15.6563 mL	31.3126 mL
5 mM	0.6263 mL	3.1313 mL	6.2625 mL
10 mM	0.3131 mL	1.5656 mL	3.1313 mL
50 mM	0.0626 mL	0.3131 mL	0.6263 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481