

## Compound T128082

## Chemical Properties

CAS No. :

Formula: C<sub>68</sub>H<sub>122</sub>O<sub>31</sub>

Molecular Weight: 1435.69

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	5-[(5-[[3-hydroxy-5-((3-hydroxy-5-[[3,4,5-trihydroxy-6-(hydroxymethyl)oxan-2-yl]oxy]decanoyl)oxy]decanoyl]oxy)decanoyl]oxy]-3-[[3,4,5-trihydroxy-6-(hydroxymethyl)oxan-2-yl]oxy]decanoyl]oxy]-3-[[3,4,5-trihydroxy-6-(hydroxymethyl)oxan-2-yl]oxy]decanoic acid is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.6965 mL	3.4826 mL	6.9653 mL
5 mM	0.1393 mL	0.6965 mL	1.3931 mL
10 mM	0.0697 mL	0.3483 mL	0.6965 mL
50 mM	0.0139 mL	0.0697 mL	0.1393 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E\_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481