

B 4148

Chemical Properties

CAS No. : 103412-37-9


Formula: C62H95N19O14S2

Molecular Weight: 1394.67

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	B 4148 is a bradykinin inhibitor.
Targets(IC50)	Others,Bradykinin Receptor

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.717 mL	3.5851 mL	7.1702 mL
5 mM	0.1434 mL	0.717 mL	1.434 mL
10 mM	0.0717 mL	0.3585 mL	0.717 mL
50 mM	0.0143 mL	0.0717 mL	0.1434 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481