


NVP-QAB205

Chemical Properties

CAS No. : 325165-07-9
 Formula: C₁₈H₂₁N₇O.xHCl
 Molecular Weight: 387.87
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	NVP-QAB205 is a Syk inhibitor.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.5782 mL	12.8909 mL	25.7818 mL
5 mM	0.5156 mL	2.5782 mL	5.1564 mL
10 mM	0.2578 mL	1.2891 mL	2.5782 mL
50 mM	0.0516 mL	0.2578 mL	0.5156 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481