


## Pyrazolo[1,5-a]pyrimidine

## Chemical Properties

CAS No. :	274-71-5
Formula:	C <sub>6</sub> H <sub>5</sub> N <sub>3</sub>
Molecular Weight:	119.12
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Pyrazolo[1,5-a]pyrimidine, with CAS No. 274-71-5, is a fragment molecule that serves as an important scaffold for molecular linking, expansion, and modification. Pyrazolo[1,5-a]pyrimidine provides a structural basis and research tool for the design and screening of novel drug candidates, and is commonly used in drug discovery, drug synthesis, and related research.
-------------	--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	8.3949 mL	41.9745 mL	83.949 mL
5 mM	1.679 mL	8.3949 mL	16.7898 mL
10 mM	0.8395 mL	4.1974 mL	8.3949 mL
50 mM	0.1679 mL	0.8395 mL	1.679 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481