

m-PEG36-OH

## Chemical Properties

CAS No. :

Formula: C73H148O37

Molecular Weight: 1617.93

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	m-PEG36-OH is a PEG-based linker for PROTACs that connects two essential ligands, facilitating the formation of PROTAC molecules and enabling selective protein degradation via the ubiquitin-proteasome system within cells.
Targets(IC50)	Others,PROTAC Linker
In vitro	PROTACs, composed of two ligands linked by a connector—one for an E3 ubiquitin ligase and the other for the target protein—utilize the intracellular ubiquitin-proteasome system to selectively degrade target proteins[1].

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.6181 mL	3.0904 mL	6.1807 mL
5 mM	0.1236 mL	0.6181 mL	1.2361 mL
10 mM	0.0618 mL	0.309 mL	0.6181 mL
50 mM	0.0124 mL	0.0618 mL	0.1236 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

## Reference

An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481