


Asperlactone

Chemical Properties

CAS No. :	76375-62-7
Formula:	C ₉ H ₁₂ O ₄
Molecular Weight:	184.191
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Asperlactone is a nematicidal, insecticidal, antibacterial, and antifungal polyketide metabolite produced from <i>A. westerdijkiae</i> .
Targets(IC50)	Others,Endogenous Metabolite

Solubility Information

Solubility	Methanol: Soluble DMSO: Soluble (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.4292 mL	27.1459 mL	54.2918 mL
5 mM	1.0858 mL	5.4292 mL	10.8584 mL
10 mM	0.5429 mL	2.7146 mL	5.4292 mL
50 mM	0.1086 mL	0.5429 mL	1.0858 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481