


Arnebin 1

Chemical Properties

CAS No. :	5162-01-6
Formula:	C ₂₁ H ₂₂ O ₆
Molecular Weight:	370.40
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	(Rac)-Arnebin 1 (beta, beta-dimethylacrylshikonin) has anti-tumor, anti-inflammatory, anti-immune deficiency and protecting liver.
Targets(IC50)	Others,HIF/HIF Prolyl-Hydroxylase

Solubility Information

Solubility	DMSO: 12.5 mg/mL (33.75 mM),Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.6998 mL	13.4989 mL	26.9978 mL
5 mM	0.540 mL	2.6998 mL	5.3996 mL
10 mM	0.270 mL	1.3499 mL	2.6998 mL
50 mM	0.054 mL	0.270 mL	0.540 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481