

EPZ005687N

Chemical Properties


CAS No. :

Formula: C₃₂H₃₇N₅O₃

Molecular Weight: 539.68

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	EPZ005687N is an inactive isomer of EPZ005687. It can also be used as the negative control for EPZ005687.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.8529 mL	9.2647 mL	18.5295 mL
5 mM	0.3706 mL	1.8529 mL	3.7059 mL
10 mM	0.1853 mL	0.9265 mL	1.8529 mL
50 mM	0.0371 mL	0.1853 mL	0.3706 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481