


Z 124-73

## Chemical Properties

CAS No. : 61930-46-9  
 Formula: C<sub>27</sub>H<sub>28</sub>N<sub>4</sub>O<sub>2</sub>S  
 Molecular Weight: 472.6  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Z 124-73 is an agent of antiviral.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.116 mL	10.5798 mL	21.1595 mL
5 mM	0.4232 mL	2.116 mL	4.2319 mL
10 mM	0.2116 mL	1.058 mL	2.116 mL
50 mM	0.0423 mL	0.2116 mL	0.4232 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

## Reference

Tonew M. Mechanism of action of antiviral hetarylhydrazones in vitro. Acta Virol. 1981 Mar;25(2):108-14. PubMed PMID: 6166184.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481