

E3 ligase Ligand 13


Chemical Properties

CAS No. :

Formula: C₃₁H₃₄N₂O₆

Molecular Weight: 530.61

Storage: Keep away from direct sunlight
 Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	E3 ligase Ligand 13 is an E3 ubiquitin ligase ligand that can be employed to create PROTACs, which are connected to the ligand for the protein through a linker. PROTACs serve as inducers of ubiquitination-mediated degradation of cancer-promoting proteins [1].
Targets(IC50)	Others,Ligands for E3 Ligase,PROTAC Linker

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.8846 mL	9.4231 mL	18.8462 mL
5 mM	0.3769 mL	1.8846 mL	3.7692 mL
10 mM	0.1885 mL	0.9423 mL	1.8846 mL
50 mM	0.0377 mL	0.1885 mL	0.3769 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Shibata N, et al. Development of Protein Degradation Inducers of Androgen Receptor by Conjugation of Androgen Receptor Ligands and Inhibitor of Apoptosis Protein Ligands. J Med Chem. 2018 Jan 25;61(2):543-575.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481