

## Exendin-4 peptide derivative acetate

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from moisture**  
Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
Actual storage temperature shall be subject to the COA.

### Biological Description

Description

Exendin-4 peptide derivative acetate is an acetate salt of Exendin-4 peptide derivative.

### Reference

US20150315260 A1

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481