

DSPE-PEG46-NH2

Chemical Properties

CAS No. :

Formula: C134H267N2O55P

Molecular Weight: 2817.51

Storage: Keep away from direct sunlight
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	DSPE-PEG46-NH2 is a PEG-based linker for PROTACs that joins two essential ligands, facilitating the formation of PROTAC molecules and enabling selective protein degradation through the ubiquitin-proteasome system within cells.
Targets(IC50)	Others,PROTAC Linker
In vitro	PROTACs consist of two ligands joined by a linker: one targets an E3 ubiquitin ligase and the other binds the target protein. They utilize the intracellular ubiquitin-proteasome system to selectively degrade target proteins [1].

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.3549 mL	1.7746 mL	3.5492 mL
5 mM	0.071 mL	0.3549 mL	0.7098 mL
10 mM	0.0355 mL	0.1775 mL	0.3549 mL
50 mM	0.0071 mL	0.0355 mL	0.071 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481