


Ylangenyl acetate

Chemical Properties

CAS No. : 90039-63-7
 Formula: C₁₇H₂₆O₂
 Molecular Weight: 262.39
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	Ylangenyl acetate is a natural product of Schisandra, Schisandraceae. The catalog number is TN5360 and the CAS number is 90039-63-7. Ylangenyl acetate can be used as a reference standard.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.8111 mL	19.0556 mL	38.1112 mL
5 mM	0.7622 mL	3.8111 mL	7.6222 mL
10 mM	0.3811 mL	1.9056 mL	3.8111 mL
50 mM	0.0762 mL	0.3811 mL	0.7622 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481