


8(14),15-Isopimaradiene-3,18-diol

Chemical Properties

| | |
|-------------------|---|
| CAS No. : | 59219-64-6 |
| Formula: | C ₂₀ H ₃₂ O ₂ |
| Molecular Weight: | 304.47 |
| Storage: | Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA. |



Biological Description

| | |
|---------------|--|
| Description | 8(14),15-Isopimaradiene-3,18-diol is a natural product for research related to life sciences. The catalog number is TN3261 and the CAS number is 59219-64-6. |
| Targets(IC50) | Others |

Preparing Stock Solutions

| | 1mg | 5mg | 10mg |
|-------|-----------|-----------|-----------|
| 1 mM | 3.2844 mL | 16.422 mL | 32.844 mL |
| 5 mM | 0.6569 mL | 3.2844 mL | 6.5688 mL |
| 10 mM | 0.3284 mL | 1.6422 mL | 3.2844 mL |
| 50 mM | 0.0657 mL | 0.3284 mL | 0.6569 mL |

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481