


Saucerneol E

Chemical Properties

CAS No. : 630096-18-3
 Formula: C₃₀H₃₄O₈
 Molecular Weight: 522.594
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	Saucerneol E is a useful organic compound for research related to life sciences. The catalog number is T124511 and the CAS number is 630096-18-3.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9135 mL	9.5677 mL	19.1355 mL
5 mM	0.3827 mL	1.9135 mL	3.8271 mL
10 mM	0.1914 mL	0.9568 mL	1.9135 mL
50 mM	0.0383 mL	0.1914 mL	0.3827 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481