


## 1-Benzyl-1,2,3,4-tetrahydroisoquinoline, sulfuric acid

## Chemical Properties

CAS No. : 1423034-30-3  
 Formula: C<sub>16</sub>H<sub>19</sub>NO<sub>4</sub>S  
 Molecular Weight: 321.4  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	1-Benzyl-1,2,3,4-tetrahydroisoquinoline, sulfuric acid, with CAS No. 1423034-30-3, is a fragment molecule that serves as an important scaffold for molecular linking, expansion, and modification. 1-Benzyl-1,2,3,4-tetrahydroisoquinoline, sulfuric acid provides a structural basis and research tool for the design and screening of novel drug candidates, and is commonly used in drug discovery, drug synthesis, and related research.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.1114 mL	15.5569 mL	31.1139 mL
5 mM	0.6223 mL	3.1114 mL	6.2228 mL
10 mM	0.3111 mL	1.5557 mL	3.1114 mL
50 mM	0.0622 mL	0.3111 mL	0.6223 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481