


AS2553627

Chemical Properties

CAS No. : 1251906-10-1
 Formula: C₁₈H₁₉N₅O
 Molecular Weight: 321.38
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	AS2553627 is a JAK inhibitor. AS2553627 prevents chronic rejection in rat cardiac allografts.
Targets(IC50)	Others,JAK

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.1116 mL	15.5579 mL	31.1158 mL
5 mM	0.6223 mL	3.1116 mL	6.2232 mL
10 mM	0.3112 mL	1.5558 mL	3.1116 mL
50 mM	0.0622 mL	0.3112 mL	0.6223 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Nakamura K, Inami M, Morio H, Okuma K, Ito M, Noto T, Shirakami S, Hirose J, Morokata T. AS2553627, a novel JAK inhibitor, prevents chronic rejection in rat cardiac allografts. Eur J Pharmacol. 2017 Feb 5;796:69-75. doi: 10.1016/j.ejphar.2016.12.025. Epub 2016 Dec 16. PubMed PMID: 27993641.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481