


## IS-741 calcium

## Chemical Properties

CAS No. :	141284-77-7
Formula:	C30H38CaF6N6O6S2
Molecular Weight:	796.86
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	IS-741 calcium is a phospholipase A2 inhibitor.
Targets(IC50)	Others,Integrin,Phospholipase

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.2549 mL	6.2746 mL	12.5493 mL
5 mM	0.251 mL	1.2549 mL	2.5099 mL
10 mM	0.1255 mL	0.6275 mL	1.2549 mL
50 mM	0.0251 mL	0.1255 mL	0.251 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481