


Hexyl hexanoate

Chemical Properties

CAS No. :	6378-65-0
Formula:	C ₁₂ H ₂₄ O ₂
Molecular Weight:	200.32
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Hexyl hexanoate is found in alcoholic beverages. Hexyl hexanoate is used in fruit flavouring. Hexyl hexanoate is present in many fruits, Parmesan cheese, alcoholic beverages and black tea. Hexyl hexanoate is a volatile component from fruit ripening.
Targets(IC50)	Others

Solubility Information

Solubility	DMSO: 25 mg/mL (124.8 mM), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.992 mL	24.9601 mL	49.9201 mL
5 mM	0.9984 mL	4.992 mL	9.984 mL
10 mM	0.4992 mL	2.496 mL	4.992 mL
50 mM	0.0998 mL	0.4992 mL	0.9984 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481