

5-Hydroxy-1,7-bis(4-hydroxyphenyl)heptan-3-yl acetate

Chemical Properties


CAS No. : 1269839-24-8

Formula: C₂₁H₂₆O₅

Molecular Weight: 358.43

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	5-Hydroxy-1,7-bis(4-hydroxyphenyl)heptan-3-yl acetate is a natural product for research related to life sciences. The catalog number is TN3115 and the CAS number is 1269839-24-8.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.7899 mL	13.9497 mL	27.8995 mL
5 mM	0.558 mL	2.7899 mL	5.5799 mL
10 mM	0.279 mL	1.395 mL	2.7899 mL
50 mM	0.0558 mL	0.279 mL	0.558 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481