


(E)-Aldosecologanin

Chemical Properties

CAS No. : 471271-55-3
 Formula: C₃₄H₄₆O₁₉
 Molecular Weight: 758.72
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	(E)-Aldosecologanin is a natural product from <i>Lonicera japonica</i> .
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.318 mL	6.590 mL	13.1801 mL
5 mM	0.2636 mL	1.318 mL	2.636 mL
10 mM	0.1318 mL	0.659 mL	1.318 mL
50 mM	0.0264 mL	0.1318 mL	0.2636 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481