


(E)-Aldosecologanin

Chemical Properties

CAS No. : 471271-55-3
 Formula: C₃₄H₄₆O₁₉
 Molecular Weight: 758.72
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

| | |
|---------------|--|
| Description | (E)-Aldosecologanin is a natural product from <i>Lonicera japonica</i> . |
| Targets(IC50) | Others |

Preparing Stock Solutions

| | 1mg | 5mg | 10mg |
|-------|-----------|-----------|------------|
| 1 mM | 1.318 mL | 6.590 mL | 13.1801 mL |
| 5 mM | 0.2636 mL | 1.318 mL | 2.636 mL |
| 10 mM | 0.1318 mL | 0.659 mL | 1.318 mL |
| 50 mM | 0.0264 mL | 0.1318 mL | 0.2636 mL |

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481