


Ko-3290

Chemical Properties

CAS No. : 79848-61-6
 Formula: C₁₉H₂₅N₃O₃
 Molecular Weight: 343.42
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	Ko-3290 is a cardioselective antagonist of β -adrenoceptors known for its antilipolytic effects in animals.
Targets(IC50)	Others, Adrenergic Receptor

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.9119 mL	14.5594 mL	29.1189 mL
5 mM	0.5824 mL	2.9119 mL	5.8238 mL
10 mM	0.2912 mL	1.4559 mL	2.9119 mL
50 mM	0.0582 mL	0.2912 mL	0.5824 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

T.H.PRINGLE, et al. The assessment of the-3-adrenoceptor blocking activity and cardioselectivity of Koe3290 in normal subjects. British Journal of Clinical Pharmacology (1987), 24(5), 581-9

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481