


(2E)-3-(4-methoxyphenyl)-1-[2,4,6-trimethoxy-3-(3-methylbut-2-en-1-yl)phenyl]prop-2-en-1-one

Chemical Properties

CAS No. :

Formula: C₂₄H₂₈O₅

Molecular Weight: 396.48

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.**Biological Description**

Description	(2E)-3-(4-methoxyphenyl)-1-[2,4,6-trimethoxy-3-(3-methylbut-2-en-1-yl)phenyl]prop-2-en-1-one is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.5222 mL	12.611 mL	25.222 mL
5 mM	0.5044 mL	2.5222 mL	5.0444 mL
10 mM	0.2522 mL	1.2611 mL	2.5222 mL
50 mM	0.0504 mL	0.2522 mL	0.5044 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481