


## 1,4-Diazepan-1-yl(2-thienyl)methanone

## Chemical Properties

CAS No. :	683274-51-3
Formula:	C10H14N2OS
Molecular Weight:	210.3
Storage:	Pure form: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	1,4-Diazepan-1-yl(2-thienyl)methanone, with CAS No. 683274-51-3, is a fragment molecule that serves as an important scaffold for molecular linking, expansion, and modification. 1,4-Diazepan-1-yl(2-thienyl)methanone provides a structural basis and research tool for the design and screening of novel drug candidates, and is commonly used in drug discovery, drug synthesis, and related research.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.7551 mL	23.7756 mL	47.5511 mL
5 mM	0.951 mL	4.7551 mL	9.5102 mL
10 mM	0.4755 mL	2.3776 mL	4.7551 mL
50 mM	0.0951 mL	0.4755 mL	0.951 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481