

RSK4-IN-1

Chemical Properties


CAS No. : 2755819-10-2

Formula: C₁₉H₂₀F₂N₄O₃

Molecular Weight: 390.38

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	RSK4-IN-1 is a compound exhibiting potent inhibition of RSK4, demonstrated by an IC ₅₀ value of 9.5 nM.
Targets(IC ₅₀)	Others,S6 Kinase

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.5616 mL	12.808 mL	25.6161 mL
5 mM	0.5123 mL	2.5616 mL	5.1232 mL
10 mM	0.2562 mL	1.2808 mL	2.5616 mL
50 mM	0.0512 mL	0.2562 mL	0.5123 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481