


CID-1792579

Chemical Properties

CAS No. : 791795-02-3
 Formula: C₂₆H₂₅N₃O₃S₂
 Molecular Weight: 491.62
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	CID-1792579 is a selective agonist of GPR55.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0341 mL	10.1705 mL	20.3409 mL
5 mM	0.4068 mL	2.0341 mL	4.0682 mL
10 mM	0.2034 mL	1.017 mL	2.0341 mL
50 mM	0.0407 mL	0.2034 mL	0.4068 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481