

(S,R,S)-AHPC-propargyl

Chemical Properties

CAS No. : 2098799-78-9


Formula: C₂₇H₃₄N₄O₅S

Molecular Weight: 526.65

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	(S,R,S)-AHPC-propargyl, a VHL ligand, is utilized in click reactions for PROTACs.
Targets(IC50)	Others,Ligands for E3 Ligase

Solubility Information

Solubility	DMSO: 100 mg/mL (189.88 mM),Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
In vivo Formulation	10% DMSO+40% PEG300+5% Tween-80+45% Saline: 3.3 mg/mL (6.27 mM),Sonication is recommended. <i>Please add the solvents sequentially, clarifying the solution as much as possible before adding the next one. Dissolve by heating and/or sonication if necessary. Working solution is recommended to be prepared and used immediately. The formulation provided above is for reference purposes only. In vivo formulations may vary and should be modified based on specific experimental conditions.</i>

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.8988 mL	9.494 mL	18.9879 mL
5 mM	0.3798 mL	1.8988 mL	3.7976 mL
10 mM	0.1899 mL	0.9494 mL	1.8988 mL
50 mM	0.038 mL	0.1899 mL	0.3798 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Wurz RP, et al. A "Click Chemistry Platform" for the Rapid Synthesis of Bispecific Molecules for Inducing Protein Degradation. J Med Chem. 2018 Jan 25;61(2):453-461.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481