


ON012380

Chemical Properties

CAS No. : 592543-24-3
 Formula: C₂₂H₂₇N₀S₀
 Molecular Weight: 465.52
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	ON012380 is a non-ATP-competitive Bcr-Abl inhibitor, potently inhibiting imatinib-resistant Bcr-Abl mutants such as T315I.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1481 mL	10.7407 mL	21.4814 mL
5 mM	0.4296 mL	2.1481 mL	4.2963 mL
10 mM	0.2148 mL	1.0741 mL	2.1481 mL
50 mM	0.043 mL	0.2148 mL	0.4296 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481