

## Dehydrojuncuenin A

## Chemical Properties


CAS No. : 1161681-26-0

Formula: C<sub>18</sub>H<sub>16</sub>O

Molecular Weight: 248.32

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

## Description

Dehydrojuncuenin A is a natural product for research related to life sciences. The catalog number is TN5694 and the CAS number is 1161681-26-0.

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.0271 mL	20.1353 mL	40.2706 mL
5 mM	0.8054 mL	4.0271 mL	8.0541 mL
10 mM	0.4027 mL	2.0135 mL	4.0271 mL
50 mM	0.0805 mL	0.4027 mL	0.8054 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481