


Tilbroquinol

Chemical Properties

CAS No. :	7175-09-9
Formula:	C ₁₀ H ₈ BrNO
Molecular Weight:	238.08
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Tilbroquinol is an antiprotozoal agent effective against amoebiasis.
Targets(IC50)	Others,Antibacterial,Parasite

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.2003 mL	21.0013 mL	42.0027 mL
5 mM	0.8401 mL	4.2003 mL	8.4005 mL
10 mM	0.420 mL	2.1001 mL	4.2003 mL
50 mM	0.084 mL	0.420 mL	0.8401 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481