


Lupanine perchlorate

Chemical Properties

CAS No. : 7400-11-5
 Formula: C₁₅H₂₅ClN₂O₅
 Molecular Weight: 348.82
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	Lupanine perchlorate is a bioactive chemical.
Targets(IC50)	Others,AChR

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.8668 mL	14.334 mL	28.6681 mL
5 mM	0.5734 mL	2.8668 mL	5.7336 mL
10 mM	0.2867 mL	1.4334 mL	2.8668 mL
50 mM	0.0573 mL	0.2867 mL	0.5734 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481