


CHK1 inhibitor

Chemical Properties

CAS No. :	2097938-64-0
Formula:	C17H21BrN4O
Molecular Weight:	377.28
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	CHK1 inhibitor (GDC-0575 analog) functions as a CHK1 inhibitor.
Targets(IC50)	Chk

Solubility Information

Solubility	DMSO: 150 mg/mL (397.58 mM),Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.6506 mL	13.2528 mL	26.5055 mL
5 mM	0.5301 mL	2.6506 mL	5.3011 mL
10 mM	0.2651 mL	1.3253 mL	2.6506 mL
50 mM	0.053 mL	0.2651 mL	0.5301 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481