


Bismaleimide

Chemical Properties

CAS No. :	13676-54-5
Formula:	C ₂₁ H ₁₄ N ₂ O ₄
Molecular Weight:	358.35
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Bismaleimide (XU 292A) has heat resistance, oxidation resistance and radiation resistance, and can be used to prepare electrical insulation materials, wear-resistant materials, reinforced plastic additives, grinding wheel adhesives, etc.
Targets(IC50)	Others

Solubility Information

Solubility	DMSO: 10 mg/mL (27.91 mM), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.7906 mL	13.9528 mL	27.9057 mL
5 mM	0.5581 mL	2.7906 mL	5.5811 mL
10 mM	0.2791 mL	1.3953 mL	2.7906 mL
50 mM	0.0558 mL	0.2791 mL	0.5581 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

- Guaresti O, García-Astrain C, Aguirresarobe RH, Eceiza A, Gabilondo N. Synthesis of stimuli-responsive chitosan-based hydrogels by Diels-Alder cross-linking `click` reaction as potential carriers for drug administration. *Carbohydr Polym.* 2018 Mar 1;183:278-286. doi: 10.1016/j.carbpol.2017.12.034. Epub 2017 Dec 16. PubMed PMID: 29352885.
- Polgar LM, Hagting E, Raffa P, Mauri M, Simonutti R, Picchioni F, van Duin M. Effect of Rubber Polarity on Cluster Formation in Rubbers Cross-Linked with Diels-Alder Chemistry. *Macromolecules.* 2017 Nov 28;50(22):8955-8964. doi: 10.1021/acs.macromol.7b01541. Epub 2017 Nov 7. PubMed PMID: 29213149; PubMed Central PMCID: PMC5707623.
- Banerjee SL, Singha NK. A new class of dual responsive self-healable hydrogels based on a core crosslinked ionic block copolymer micelle prepared via RAFT polymerization and Diels-Alder "click" chemistry. *Soft Matter.* 2017 Dec 6;13(47):9024-9035. doi: 10.1039/c7sm01906h. PubMed PMID: 29177283.
- Berto P, Grelier S, Peruch F. Telechelic Polybutadienes or Polyisoprenes Precursors for Recyclable Elastomeric Networks. *Macromol Rapid Commun.* 2017 Nov;38(22). doi: 10.1002/marc.201700475. Epub 2017 Oct 12. PubMed PMID: 29024134.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481