


20(29)-Lupene-3,23-diol

Chemical Properties

CAS No. :	32451-85-7
Formula:	C ₃₀ H ₅₀ O ₂
Molecular Weight:	442.72
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	20(29)-Lupene-3,23-diol is a natural product for research related to life sciences. The catalog number is TN2796 and the CAS number is 32451-85-7.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2588 mL	11.2938 mL	22.5876 mL
5 mM	0.4518 mL	2.2588 mL	4.5175 mL
10 mM	0.2259 mL	1.1294 mL	2.2588 mL
50 mM	0.0452 mL	0.2259 mL	0.4518 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481