


Geranyl Tiglate

Chemical Properties

CAS No. :	8000-46-2
Formula:	C ₁₅ H ₂₄ O ₂
Molecular Weight:	236.36
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Geranyl Tiglate (tiglic acid geranyl ester) is one of the most important natural perfume.
Targets(IC50)	Estrogen Receptor/ERR,Others

Solubility Information

Solubility	DMSO: 45 mg/mL (190.39 mM),Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.2308 mL	21.1542 mL	42.3083 mL
5 mM	0.8462 mL	4.2308 mL	8.4617 mL
10 mM	0.4231 mL	2.1154 mL	4.2308 mL
50 mM	0.0846 mL	0.4231 mL	0.8462 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Sienkiewicz M, et al. The antibacterial activity of geranium oil against Gram-negative bacteria isolated from difficult-to-heal wounds. Burns. 2014 Aug;40(5):1046-51.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481