


Torilolone

Chemical Properties

CAS No. :	13018-09-2
Formula:	C ₁₅ H ₂₄ O ₃
Molecular Weight:	252.35
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Torilin and torilolone show hepatoprotective effects on tacrine-induced cytotoxicity in human liver-derived Hep G2 cells, the EC ₅₀ values are 20.6 +/- 1.86 and 3.6 +/- 0.1 microM, respectively; and silybin as a positive control shows an EC ₅₀ value of 69.0 +/- 3.4 microM.
Targets(IC ₅₀)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.9628 mL	19.8138 mL	39.6275 mL
5 mM	0.7926 mL	3.9628 mL	7.9255 mL
10 mM	0.3963 mL	1.9814 mL	3.9628 mL
50 mM	0.0793 mL	0.3963 mL	0.7926 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481