


Seridopidine

Chemical Properties

CAS No. :	883631-51-4
Formula:	C ₁₄ H ₂₀ FNO ₂ S
Molecular Weight:	285.38
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Seridopidine (ACR343) is a modulator of dopaminergic activity that is used as an oral therapy for schizophrenia, Parkinson's disease, and Tourette's syndrome.
Targets(IC ₅₀)	Dopamine Receptor

Solubility Information

Solubility	DMSO: 45 mg/mL (157.68 mM), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.5041 mL	17.5205 mL	35.041 mL
5 mM	0.7008 mL	3.5041 mL	7.0082 mL
10 mM	0.3504 mL	1.752 mL	3.5041 mL
50 mM	0.0701 mL	0.3504 mL	0.7008 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Sonesson, et al. 3,5-Disubstituted phenyl-piperidines as modulators of dopamine neurotransmission, their preparation, pharmaceutical compositions, and use in therapy. WO2007042295A1.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481