

(2S,3R)-LP99

## Chemical Properties


CAS No. : 1808948-28-8

Formula: C<sub>26</sub>H<sub>30</sub>ClN<sub>3</sub>O<sub>4</sub>S

Molecular Weight: 516.05

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	(2S,3R)-LP99 is the less active enantiomer of LP99.
Targets(IC50)	Others,Epigenetic Reader Domain

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9378 mL	9.689 mL	19.378 mL
5 mM	0.3876 mL	1.9378 mL	3.8756 mL
10 mM	0.1938 mL	0.9689 mL	1.9378 mL
50 mM	0.0388 mL	0.1938 mL	0.3876 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481