

Cathepsin C-IN-3

Chemical Properties

CAS No. :

Formula: C₂₈H₂₁F₃N₆O₅

Molecular Weight: 546.57

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Cathepsin C-IN-3 is a potent inhibitor of histone C (IC ₅₀ : 61.79 nM) and also inhibits THP-1 cells (IC ₅₀ : 101.5 nM) and U937 cells (IC ₅₀ : 86.5 nM).
Targets(IC ₅₀)	Others,Cysteine Protease

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.8296 mL	9.148 mL	18.2959 mL
5 mM	0.3659 mL	1.8296 mL	3.6592 mL
10 mM	0.183 mL	0.9148 mL	1.8296 mL
50 mM	0.0366 mL	0.183 mL	0.3659 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481