

Compound T130675

Chemical Properties


CAS No. :

Formula: C₁₈H₂₀N₂O₉

Molecular Weight: 408.36

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	(2R,3S,4R)-5-(2,4-dioxo-1,2,3,4-tetrahydropyrimidin-1-yl)-4-hydroxy-2-(hydroxymethyl)oxolan-3-yl 4-hydroxy-2-methoxy-6-methylbenzoate is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.4488 mL	12.2441 mL	24.4882 mL
5 mM	0.4898 mL	2.4488 mL	4.8976 mL
10 mM	0.2449 mL	1.2244 mL	2.4488 mL
50 mM	0.049 mL	0.2449 mL	0.4898 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481