


Stachartin E

Chemical Properties

CAS No. :	1978388-58-7
Formula:	C ₃₃ H ₄₁ N ₂ O ₆
Molecular Weight:	547.68
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Stachartin E is a natural product from the cultures of the tin mine tailings-associated fungus <i>Stachybotrys chartarum</i> .
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.8259 mL	9.1294 mL	18.2588 mL
5 mM	0.3652 mL	1.8259 mL	3.6518 mL
10 mM	0.1826 mL	0.9129 mL	1.8259 mL
50 mM	0.0365 mL	0.1826 mL	0.3652 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481