


## 4-Aminophenylacetic acid

## Chemical Properties

CAS No. :	1197-55-3
Formula:	C <sub>8</sub> H <sub>9</sub> NO <sub>2</sub>
Molecular Weight:	151.16
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	4-Aminophenylacetic acid ,with CAS No. 1197-55-3, is a fragment molecule that serves as an important scaffold for molecular linking, expansion, and modification. 4-Aminophenylacetic acid provides a structural basis and research tool for the design and screening of novel drug candidates, and is commonly used in drug discovery, drug synthesis, and related research.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	6.6155 mL	33.0775 mL	66.1551 mL
5 mM	1.3231 mL	6.6155 mL	13.231 mL
10 mM	0.6616 mL	3.3078 mL	6.6155 mL
50 mM	0.1323 mL	0.6616 mL	1.3231 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481