


PERK-IN-3

Chemical Properties

CAS No. :	1337532-08-7
Formula:	C ₂₂ H ₁₆ F ₂ N ₄ O ₂
Molecular Weight:	406.38
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	PERK-IN-3 is a potent inhibitor of PERK (IC ₅₀ of 7.4 nM).
Targets(IC ₅₀)	ERK,PERK

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.4608 mL	12.3038 mL	24.6075 mL
5 mM	0.4922 mL	2.4608 mL	4.9215 mL
10 mM	0.2461 mL	1.2304 mL	2.4608 mL
50 mM	0.0492 mL	0.2461 mL	0.4922 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Axten JM, et al. Discovery of 7-methyl-5-(1-[[3-(trifluoromethyl)phenyl]acetyl]-2,3-dihydro-1H-indol-5-yl)-7H-pyrrolo[2,3-d]pyrimidin-4-amine (GSK2606414), a potent and selective first-in-class inhibitor of protein kinase R (PKR)-like endoplasmic reticulum kinase (PERK). J Med Chem. 2012 Aug 23;55(16):7193-207.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481