

BRD73842

## Chemical Properties


CAS No. : 1643982-19-7

Formula: C<sub>28</sub>H<sub>33</sub>FN<sub>4</sub>O<sub>4</sub>

Molecular Weight: 508.58

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	BRD73842 is an inhibitor of Plasmodium PI4K.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9663 mL	9.8313 mL	19.6626 mL
5 mM	0.3933 mL	1.9663 mL	3.9325 mL
10 mM	0.1966 mL	0.9831 mL	1.9663 mL
50 mM	0.0393 mL	0.1966 mL	0.3933 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481